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EN | INSTRUCTIONS FOR USE

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ZES-CA-D-5117-C-Rev-01
Enero 2026

ZIADERMA

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ZIADERMA

1. Descripción

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ZIADERMA

1. Description

ZIADERMA is a sterile medical device that is injectable, apyrogenic, and reabsorbable. It contains reticulated hyaluronic acid from bacterial fermentation. ZIADERMA comes in the form of a clear, colourless, consistent gel, contained in pre-filled disposable syringes with Luer-lock connector.

2. Composition

- ZIADERMA SOFT: Concentrated reticulated hyaluronic acid. Sterile saline solution: q.s.
- ZIADERMA MEDIUM: Medium concentration reticulated hyaluronic acid. Sterile saline solution: q.s.
- ZIADERMA DEEP: High concentration reticulated hyaluronic acid. Sterile saline solution: q.s.

3. Action mechanism

ZIADERMA's action consists of boosting the volume of the dermal-epidermal tissue, based on the natural capacity of the hydrophilic molecules of hyaluronic acid to bond with an amount of water many times greater than their own weight. This makes it possible to fill the intradermal spaces and integrate the intercellular matrix, firming up the tissue. The reticulation of hyaluronic acid inside the product makes the filling effect more stable and longer lasting.

4. Instructions for use

ZIADERMA is intended to be a temporary filler to correct skin defects such as wrinkles or scars. It is also suitable for use on pathological conditions such as lipoatrophy linked to HIV infection and anti-retroviral treatments. Particularly, the following use is recommended:

- ZIADERMA SOFT to correct lines and surface wrinkles, from fine to medium lines. Defects and particularly indicated areas: laughter lines, bar code lip lines, wrinkles around the eyes, scarring from acne, chicken pox.
- ZIADERMA MEDIUM to correct medium-surface wrinkles on the face. Defects and particularly indicated areas: nasolabial wrinkles, nose, lips, frown lines.
- ZIADERMA DEEP to correct deep facial wrinkles and furrows. Defects and particularly indicated areas: nasolabial furrows, marionette mouth wrinkles, temples, reshaping facial contours, lips.

5. Contraindications

ZIADERMA must not be used: in injection form in the periorbital area (eyelids, dark circles); by injection into blood vessels; to increase breast tissue or for insertion into bones, tendons, ligaments or muscles. ZIADERMA must not be administered to patients with: pre-existing or ongoing autoimmune diseases, or which are being treated with immunotherapy drugs; verified tendency to develop keloids or hypertrophic scars, or pigmentation disorders; serious multiple allergies; infections and inflammatory processes in or near the area to be treated; coagulation problems; recognised hypersensitivity to hyaluronic acid. ZIADERMA must not be used by pregnant or breastfeeding women or on persons under the age of 18 years.

6. Warnings

ZIADERMA is only intended for subcutaneous or intradermal injection, and must not be injected by intravascular means. When injecting, there is a potential risk of inadvertent injection of the material into a blood vessel. In rare cases, this could lead to vascular occlusion, abscess, necrosis or embolism. To ensure that the needle is not inside a blood vessel during normal use of the medical device, gently aspirate with the syringe before proceeding with the injection. ZIADERMA is made to be used on a single patient and during a single session. ZIADERMA is supplied as a sterile, single-use product and does not require re-sterilising. Not for re-use. Use ZIADERMA immediately after opening the blister pack and dispose of the syringe immediately after use, even if not all of the contents are used. Do not use the product after the expiry date shown on the packaging or if the packaging is damaged, open, or has been stored improperly. Do not use ZIADERMA together with other injectable products, or mix it with other products before use. Do not use ZIADERMA if the skin is red, swollen or shows signs of discolouration; do not use ZIADERMA when the patient is suffering from inflammation, infection or fever. Patients who take drugs to slow blood clotting may show signs of injury or pronounced bleeding at the injection site. Patients being treated with anticoagulants must suspend use of these, in agreement with their own doctor, at least 14 days before receiving injections of ZIADERMA.

INCOMPATIBILITY: there is a recognised incompatibility between hyaluronic acid and quaternary ammonium salts, such as benzalkonium chloride. Avoid all contact between ZIADERMA and these substances and do not use medical-surgical material containing quaternary ammonium salts. Keep the product out of the reach of children.

7. Precautions for use

Only use ZIADERMA according to these instructions for use. ZIADERMA is a clear, colourless gel. If the contents of the syringe appear altered or to contain particles, do not use, and inform VIDA S.r.l. Immediately at the address below. Do not use the device if the syringe or the blister in which it is packaged are damaged. Follow all normal precautions of use associated with procedures for intradermal injection. Injections are associated with an intrinsic risk of infection. The skin must be perfectly clean and disinfected before treatment. Do not use on patients with ongoing infections or inflammatory processes near or on the site of the injection. It is recommended not to use ZIADERMA in areas where other fillers have been used, and do not use it at the same time as any other injectable product. Advise patients not to use make-up in the 24 hours following injection and inform them that, during the remission period, or until all swelling has disappeared, they must avoid exposing the treated area to sunlight, UV lamps, or cold. If using laser treatment, chemical peels or dermoabrasion after the use of ZIADERMA, there is a risk of an inflammatory reaction. For the same reason, ZIADERMA must only be used when the outcomes of any previous treatment have completely ceased.

Syringes and needles are a potential biohazard after use. Handle with care and dispose of in compliance with approved medical practice.

8. Adverse reactions

Before commencing treatment, the doctor must inform the patient of the existence of potential secondary effects, linked to intradermal filling, and which can occur immediately or some time after. Following the use of ZIADERMA, it is possible for adverse reactions to occur related to the injection itself. These include: pain, itching, alterations in pigment and tissue density in the injection site. They normally go away on their own, within a week. At times, small oedema may occur, but these are reabsorbed in just a few days. Among these adverse reactions related to the use of the product, in consideration of bibliographical studies in existing literature, and the scientific knowledge of the producer at the date of issue of this informative leaflet, it is possible to highlight the following: inflammatory response, itching, pain to the touch after the injection, lasting for up to a maximum of one week; bleeding or bruising, in rare circumstances; follicular pustules-pimples that may form

ZIADERMA

up to 4 weeks after treatment, with a duration of two weeks; allergic reactions to hyaluronic acid; hardening or appearance of nodules in the injection site; in very rare circumstances, granuloma or significant inflammatory responses, hives, necrosis, rashes, which may appear with a delay of 2-4 weeks.

IMPORTANT: patients must report any serious adverse reactions or reactions not listed in this illustrative leaflet to their medical practitioner, and the practitioner, as well as prescribing suitable treatment for same, is obliged to inform distributors and the manufacturer VIDA S.r.l. without any undue delay.

9. Use

Before starting treatment, the doctor must inform the patient about the composition, characteristics, precautions for use, warnings, contraindications, and potential adverse reactions to this substance. The doctor must also proceed to assess the patient's medical history and assess any previous immunological or allergy problems with the patient, ensuring compatibility with the use of ZIADERMA. Finally, the doctor must assess the skin areas being treated to ensure patient suitability for the implant and actual usefulness of the treatment. ZIADERMA may only be administered by qualified medical personnel. For the correct end result of the treatment, it is essential to proceed with the correct injection technique. Before injection, disinfect the area to be treated, using alcohol or other antiseptic solution. This product must be injected using a sterile syringe with fine needle. For ZIADERMA SOFT we recommend a 30G needle, for ZIADERMA MEDIUM a 30G needle, for ZIADERMA DEEP a 27G or 25G needle. It is up to the doctor to decide the needle to use and the injection method. The quantity of product to be injected is established according to the depth of the wrinkle or the skin lesion, and can be easily dosed, thanks to the use of graduated syringe. Immediately after the injection, it is necessary to press firmly and massage the treated area in order to aid the even distribution of the product. Before every procedure, it is necessary to carefully assess the symmetry of the chosen facial areas, planning to inject the same amount of material in the right and left sides of the face, and in symmetrical points. Immediately before adding the filler, remove the syringe from the protective blister, unscrew the cap and fit the needle, screwing the cone onto the Luer-lock connection of the syringe, until you notice a some resistance. To use ZIADERMA safely, it is important for the needle to be suitable and correctly fitted to the syringe, since incorrect assembly may cause the needle to come out of the syringe during the injection process. Only remove the protective cap from the needle just before performing the injection. Do not use non-sterile, bent, or deformed needles or needles from damaged packaging. Before injecting, expel the air from the syringe until a droplet of product comes out on the tip of the needle. Introduce the needle at a 30-degree angle, following the exact path of the wrinkle or hollow, and inject the product into the middle dermis. If hyaluronic acid is introduced too deeply, the final effect will last for less time, while if injected too closely to the surface, it is possible for nodules or granulomas to form. Never apply too much pressure on the syringe. Inject ZIADERMA while slowly retracting the needle, depositing an even amount of product along the chosen site. Interrupt pressure on the plunger of the syringe before removing the needle from the skin to make sure no material comes out of the injection point. In case of immediate swelling, apply a cold compress for a short period. Please note: the graduation on the syringe does not serve as a measurement, but merely as a means to visually control the material injected.

10. Dosage

The dosage depends on the condition of the patient's skin and the depth of the area to be treated. It is possible to achieve excellent results with regular treatments, taking into account the different factors in deciding how long the injection remains effective, such as lifestyle, physical activity, exposure to sunlight or sun lamps, age, and suitable filling technique.

11. Packaging

A package of ZIADERMA contains: a syringe, pre-filled with the product, with 1mL graduations, for ZIADERMA SOFT, ZIADERMA MEDIUM, ZIADERMA DEEP.

Each syringe comes with Luer-Lock connection inside a blister pack, and is provided with an adhesive label. A set of two adhesive labels per syringe, with the batch number and expiry date, one to be applied to the patient's medical records and the other to be delivered to the patient. It is obligatory to follow these procedures in order to ensure product traceability. Instructions for use, to be read with care before using the device. Two 27G needles for ZIADERMA DEEP and two 30G needles for ZIADERMA SOFT and ZIADERMA MEDIUM.

12. Storage

Read the expiry date for the product on the adhesive label on the syringe or on the packet. Store in a dry place, away from direct sunlight, and at a temperature from 5°C-25°C.

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UBGEN® Srl